




#### A. POWER PACK





A.1 Carabiner  
A.2 Battery Indicator  
A.3 Connector  
A.4 Cable strap  
A.5 Cable

#### B. °M WARMER

B.1 Adhesive for fixation  
B.2 Luer to IV-fluid  
B.3 Function Indicator arrow  
B.4 Alarm Indicator  
B.5 Luer to IV-catheter  
B.6 Pinch Clamp

#### C. CHARGER

C.1 Connector to Power Pack  
C.2 Charging status  
C.3 Cable



°MEQU

Curious?

See videos and other supportive material online at [www.mequ.dk](http://www.mequ.dk)

°MEQU  
Ole Maalaes Vej 3  
2200 Copenhagen N  
Denmark  
MQ1-W032-3.0  
2018-01-20

°MEQU

°M Warmer System